

# Automated Gating of Flow Cytometry Data using the Bioconductor **openCyto** Framework

*Nichole Monhait, MPH Candidate  
Colorado School of Public Health, Colorado State University*

*2019-05-14*



# Contents

|          |   |           |
|----------|---|-----------|
| <b>1</b> | <b>What's inside?</b>                             | <b>5</b>  |
| <b>2</b> | <b>Getting Started</b>                            | <b>7</b>  |
| 2.1      | Required Packages and Installation . . . . .      | 7         |
| <b>3</b> | <b>Working with your Manual Gating Scheme</b>     | <b>9</b>  |
| 3.1      | Read in flowJo file . . . . .                     | 9         |
| 3.2      | Parse FCS files . . . . .                         | 10        |
| 3.3      | Visualize and Verify . . . . .                    | 10        |
| <b>4</b> | <b>Create .csv</b>                                | <b>13</b> |
| 4.1      | .csv Gating Template Structure . . . . .          | 13        |
| 4.2      | Creating the Template . . . . .                   | 15        |
| <b>5</b> | <b>Automate Gating</b>                            | <b>17</b> |
| 5.1      | Load .csv into R . . . . .                        | 17        |
| 5.2      | Read in raw FCS files . . . . .                   | 17        |
| 5.3      | Apply Gating . . . . .                            | 17        |
| 5.4      | Plot Automated Gating . . . . .                   | 17        |
| 5.5      | Population Statistics . . . . .                   | 19        |
| <b>6</b> | <b>Plot Customization</b>                         | <b>21</b> |
| 6.1      | Adjusting plots . . . . .                         | 21        |
| <b>7</b> | <b>Appendix</b>                                   | <b>25</b> |
| 7.1      | Function Definitions . . . . .                    | 25        |
| 7.2      | R Object Definitions . . . . .                    | 25        |
| 7.3      | CSU Flow Cytometer Information . . . . .          | 25        |
| 7.4      | Transform data for better visualization . . . . . | 25        |
| <b>8</b> | <b>References</b>                                 | <b>27</b> |



# Chapter 1

## What's inside?

Flow cytometry is a method used to gain understanding of cell samples and populations by quantifying scattered and emitted fluorescent light (7). Signals are captured and analyzed through use of software programs. Flow cytometry analysis consists of gating, a method that dictates which cells will be further analyzed and which will not (7). Current methods for flow cytometry gating involve manually drawing gates. This process is both time consuming and costly, making automated gating procedures an appealing option. The **openCyto** package allows users to take manually gated data from flowJo, reproduce those gates in R, and eventually automate the gating process. The goal of this tutorial is to take the user through the process of automated gating analysis.

This tutorial will be useful to anyone who has done manual gating on a sample and wishes to automate the same procedure on additional samples in the future.

The example data used in this tutorial is from Colorado State University's Microbiology, Immunology, and Pathology Department. For documentation on the flow cytometer used at CSU, please see the Appendix. Alternatively, you can input your own data using the filetypes described in Chapter 3.



## Chapter 2

# Getting Started

Here is an overview of the process to automate flow cytometry data using R's **openCyto** and what you will need to successfully automate your own flow cytometry analysis. The general steps to accomplish this are as follows:

1. Read in a manually gated flowJo workspace in .wsp file format.
2. Parse raw FCS files from the read in workspace.
3. Visualize the manual gating template and resulting gates to verify gating scheme.
4. Create and read in a .csv gating template.
5. Read in raw FCS files to be gated.
6. Automate gating.
7. Visualize automated gating template and gates to verify gating scheme.
8. Extract population statistics and relevant information.

This process is completed primarily with the **openCyto** package but calls upon other packages within the Bioconductor **openCyto** framework. Packages needed to complete this tutorial are listed at the end of this chapter. Descriptions of each function and R object used for this analysis can be found in the Appendix.

## 2.1 Required Packages and Installation

### 2.1.1 Package descriptions

Below is a description of each package used in this analysis. Code to install and use these packages will follow. Package descriptions were taken from Bioconductor and CRAN.

| Package Name  |   |
|---------------|---|
| openCyto      | This package is designed  |
| flowWorkspace | This package allows you to import basic flowJo workspaces into BioConductor and replicate the gating fr |
| flowCore      | Pro   |
| flowStats     | Methods and functiona   |
| flowClust     | Robu  |
| data.table    | Fast aggregation of large data (e.g. 100GB in RAM), fast ordered joins, fast add/modify/delete          |

### 2.1.2 Directory set-up

The working directory is a location on your computer where R will save and read files from. Download the folder titled “tutorial” which contains all of the content needed to complete this tutorial and place it on your desktop. To set this folder as the current working directory, use the following code. Check this by next calling, `getwd()`. The result of `getwd()` should be a pathname to the tutorial folder on your desktop. Once this is set, all of the code within this tutorial will run without needing to change pathnames.

```
setwd("/Users/monhait/Desktop/tutorial")
getwd()
```

### 2.1.3 Installation

Install the following libraries into a new R script. As you will see below, this tutorial uses the development version of `openCyto`. It is important to use the development version of `openCyto` to remain up to date on any changes made by the developers of `openCyto`. Use the following to ensure the correct packages are installed. Installation will only need to be done once.

*To install*

```
if (!requireNamespace("BiocManager", quietly = TRUE))
  install.packages("BiocManager")
BiocManager::install()

BiocManager::install(c("openCyto", "flowWorkspace", "flowCore",
                      "flowStats", "flowClust"))

install.packages("data.table")
devtools::install_github("RGLab/openCyto", ref = "trunk")
```

RStudio may also prompt you to download XQuartz and XCode based on your computer type, so it may be a good idea to go ahead and also download both.

### 2.1.4 Load packages

Although installation only needs to be done once, packages will need to be reloaded each time you open an R session. At the beginning of each session, run the following code.

*To load*

```
library(openCyto)
library(flowWorkspace)
library(data.table)
library(flowCore)
library(flowStats)
library(flowClust)
```



## Chapter 3

# Working with your Manual Gating Scheme

The first step in this process is to bring a pre-existing flowJo file into R to recreate the gating environment. This chapter will detail the following:

1. Read in flowJo *.wsp* file
2. Parse FCS files
3. Visualize and verify manual gates

### 3.1 Read in flowJo file

Within flowJo, transformation, compensation, and gating can be saved as either *.xml* or *.wsp* filetypes. This tutorial will only detail steps from a *.wsp* filetype saved from flowJo. Note that many other tutorials begin from a *.xml* filetype. Saving analysis within flowJo is detailed here.

Your *.wsp* file will contain cell samples and groups to be added to the Workspace in R. The R Workspace is where all gates, analyses, and compensation matrices applied to samples within flowJo will be replicated and can be added to other samples in R. Importantly, the *.wsp* will not contain your FCS files. Rather, the path to your files will be saved within the *.wsp* and will be referenced later in R.

Before you begin, be sure you have loaded the required packages outlined in the previous chapter.

Once all packages are loaded, save the *.wsp* filetype path as an R object. Here, it is called `wsfile`, but you can name your R objects in any way. Next, use the `openWorkspace()` function with your R object created in the prior step to open the *.wsp* file in R. Save this as a new R object of any name. Here it is named `ws`. This new R object will be of `flowJoWorkspace` class.

Here is an example of saving and opening your *.wsp* filetype in R. Please ensure that `ws` is saved as a `flowWorkspace` object containing groups of samples before proceeding. This can be verified by looking in the Environment tab in the upper right pane of your screen. Next to each saved R object will list the class that the object is saved as.

```
wsfile <- "../tutorial/group1_v_group2.wsp"
```

```
ws <- openWorkspace(wsfile)
```

```
print(ws)

## FlowJo Workspace Version 20.0
## File location: ../tutorial
## File name: group1_v_group2.wsp
## Workspace is open.
##
## Groups in Workspace
##      Name Num.Samples
## 1 All Samples      10
## 2 Samples         10
```

## 3.2 Parse FCS files

The next step is to read in raw FCS files. FCS files contain data from the flow cytometer. Standards for FCS files are listed here.

Raw FCS files are read using the `parseWorkspace()` function. This function will read the FCS files and transform, compensate, and gate according to parameters defined from the `.wsp` flowJo workspace, which is now saved as an R object of class `flowWorkspace`. The `parseWorkspace()` call requires the object that results from running `openWorkspace()`. Here, we named this object `ws`. The function `parseWorkspace()` also requires the name of the samples to read in. To list sample names, use the `getSampleGroups()` function on your `flowWorkspace` class object. You can then choose the sample names based on the results of `getSampleGroups()`.

A new R object, here named `gating_set`, is then created and will be a `GatingSet` object. The `isNcdf = TRUE` call saves this output to disk rather than into memory because the FCS files are large. The `sampNloc = 'sampleNode'` tells R where to retrieve the FCS file names. The options for `sampNloc` can be either “keyword” or “sampleNode”.

Here is an example of parsing FCS files. As this function runs, you will see several messages appear as the FCS files are loaded and the manual gating scheme is replicated. After this, `attributes()` is used to examine the data.

```
gating_set <- parseWorkspace(ws, name = "Samples",
                             path = "../tutorial/group1_v_group2",
                             isNcdf = TRUE, sampNloc = 'sampleNode')
```

```
## windows version of flowJo workspace recognized.
## version X
```

```
attributes(gating_set)
```

## 3.3 Visualize and Verify

It is helpful to now visualize both the gating template and gates on a subset of the data in order to verify the gating scheme. This will ensure consistency between the flowJo workspace and the manual gates recreated in R. First, save a subset of the `GatingSet` object type, here named `gating_set`. The following saves the first FCS file of `gating_set` as `gh`, which will be a `GatingHierarchy` object type. Since each FCS file corresponds to an individual experiment, this saves the first experiment of the group.


```
gh <- gating_set[[1]]
print(gh)
```

```
## Sample: X_group1_1
## GatingHierarchy with 51 gates
```

### 3.3.1 plot()

The `plot()` function will visualize the current gating hierarchy when applied to an object of class `GatingHierarchy`. This can be done for the entire gating hierarchy or a specific population as seen below.


```
plot(gh)
```



### 3.3.2 plotGate()

The `plotGate()` function will gate the designated subset of your data according to parameters replicated from flowJo. This must be called on an object of class `GatingHierarchy`.

```
flowWorkspace.par.set("plotGate", list(xlim = "data",
                                       ylim = "data"))
plotGate(gh, xbin = 128)
```



**\*\*Note the use of `flowWorkspace.par.set()` and `xbin` here. Chapter 5 of this tutorial will discuss these plot customizations.**

**\*\*Note the use of `flowWorkspace.par.set()` and `xbin` here. Chapter 5 of this tutorial will discuss these plot customizations.**

# Chapter 4

## Create .csv

The creation of a *.csv* gating template is arguably the most important step to automating flow cytometry analysis. The *.csv* template that you create will tell **openCyto** how to gate your data based on parameters listed within the file. Included with this tutorial is a partial *.csv* template that can be used to gate the sample data. Generally, *.csv* templates will look like this:

| alias       | pop | parent    | dims              | gating_method | gating_args                         | collapseDataForGating | groupBy | preprocessing_method | preprocessing_args |
|-------------|-----|-----------|-------------------|---------------|-------------------------------------|-----------------------|---------|----------------------|--------------------|
| nonDebris   | +   | root      | FSC-A             | mindensity    |                                     |                       |         |                      |                    |
| singlets    | +   | nonDebris | FSC-A,FSC-H       | singleGate    | prediction_level=0.99,maxit=20      |                       |         |                      |                    |
| lymph       | +   | singlets  | FSC-A,SSC-A       | flowClust     | K=2,target=q[1e5,5e4],quantile=0.99 |                       |         | prior_flowClust      |                    |
| viable      | -   | lymph     | Zombie-Nir-A      | tailgate      | tol=0.05                            |                       |         |                      |                    |
| CD3         | +   | viable    | Alexa Fluor 532-A | mindensity    |                                     |                       |         |                      |                    |
| CD4         | +   | CD3       | BV480-A           | mindensity    |                                     |                       |         |                      |                    |
| CD8         | +   | CD3       | BV570-A           | mindensity    |                                     |                       |         |                      |                    |
| CD45RB(CD8) | +   | CD8       | Pacific Blue-A    | tailgate      |                                     |                       |         |                      |                    |
| CD153(CD8)  | +   | CD8       | BV711-A           | tailgate      |                                     |                       |         |                      |                    |
| CTLA4(CD8)  | +   | CD8       | BV605-A           | tailgate      |                                     |                       |         |                      |                    |
| CXCR1(CD8)  | +   | CD8       | PE-A              | tailgate      |                                     |                       |         |                      |                    |
| CD44(CD8)   | +/- | CD8       | APC-Fire 750-A    | mindensity    |                                     |                       |         |                      |                    |

### 4.1 .csv Gating Template Structure

In the gating template, each row corresponds to a single cell population and the method used to gate that population. When read into R, the *.csv* will direct gating based on parameters listed in each row and column. The *.csv* must contain 10 predefined columns as seen here:

| alias | pop | parent | dims | gating_method | gating_args | collapseDataForGating | groupBy | preprocessing_method | preprocessing_args |
|-------|-----|--------|------|---------------|-------------|-----------------------|---------|----------------------|--------------------|
|-------|-----|--------|------|---------------|-------------|-----------------------|---------|----------------------|--------------------|

#### 4.1.1 alias

The first column must be titled **alias**. This is where you will put your cell population names. Remember, each row corresponds to a single cell population. Population names in the **alias** column must be unique.

| alias       |
|-------------|
| nonDebris   |
| singlets    |
| lymph       |
| viable      |
| CD3         |
| CD4         |
| CD8         |
| CD45RB(CD8) |
| CD153(CD8)  |
| CTLA4(CD8)  |
| CXCR1(CD8)  |
| CD44(CD8)   |

### 4.1.2 pop

The second column must be titled **pop**. This column will contain a + or - to designate which subset or quadrant will be gated. A + will gate the positive subset while a - will gate the negative. This column can only contain strings of + and -, so do not use any characters as separators for quadrant gates.

| pop |
|-----|
| +   |
| +   |
| +   |
| -   |
| +   |
| +   |
| +   |
| +   |
| +   |
| +   |
| +   |

### 4.1.3 parent

The third column must be titled **parent**. This column refers to the parent cell population, or where the current cell population originates from. Similar to the **alias** column, **parent** names must be unique. This

column cannot contain any commas, otherwise `openCyto` will assume the population has multiple parents and you will get an error message.

| parent    |
|-----------|
| root      |
| nonDebris |
| singlets  |
| lymph     |
| viable    |
| CD3       |
| CD3       |
| CD8       |
| CD8       |
| CD8       |
| CD8       |
| CD8       |
| CD8       |

#### 4.1.4 Remaining template columns

`dims`- channel or marker names for gating

`gating_method`- gating function (supported options listed above)

`-quadrantGate`

`-rangeGate`

`-quantileGate`

`-mindensity`

`-tailgate`

`-cytokine`

`-flowClust`

`-boundary`

`-singletGate`

`-transitional`

`-polyfunctionalityGate`

`-flowDensity`

`gating_args`- arguments to be passed to gating function `collapseDataforGating`- data is collapsed and replicated across all samples

`groupBy`- used to group samples into unique combinations

`preprocessing_method`- preprocessing function

`preprocessing_args`- arguments for preprocessing function

## 4.2 Creating the Template

The gating template can be created manually or assisted by the use of the `templateGen()` function. The function `templateGen()` will input the `alias`, `pop`, `parent`, and `dims` columns and the rest must be completed manually. To use `templateGen()`, you must input a `GatingHierarchy` object. In this example, that is `gh`, the subset created from `gating_set`.

```
gt <- templateGen(gh)
head(gt)
```

```
##          alias          pop          parent
## 1   Singlets   Singlets          root
## 2     Live     Live          /Singlets
## 3 Leukocytes Leukocytes    /Singlets/Live
## 4     CD3+     CD3+    /Singlets/Live/Leukocytes
## 5     CD8+     CD8+    /Singlets/Live/Leukocytes/CD3+
## 6 CXCR1-     CXCR1- /Singlets/Live/Leukocytes/CD3+/CD8+
##          dims gating_method gating_args
## 1          FSC-A,FSC-H      <NA>      <NA>
## 2      Comp-Zombie Nir-A      <NA>      <NA>
## 3          FSC-A,SSC-A      <NA>      <NA>
## 4 Comp-Alexa Fluor 532-A,SSC-H      <NA>      <NA>
## 5      Comp-BV570-A,SSC-H      <NA>      <NA>
## 6          Comp-PE-A      <NA>      <NA>
## collapseDataForGating groupBy preprocessing_method preprocessing_args
## 1          <NA>      <NA>      <NA>      <NA>
## 2          <NA>      <NA>      <NA>      <NA>
## 3          <NA>      <NA>      <NA>      <NA>
## 4          <NA>      <NA>      <NA>      <NA>
## 5          <NA>      <NA>      <NA>      <NA>
## 6          <NA>      <NA>      <NA>      <NA>
```

The auto-filled template will generate within the R Console and can then be saved locally with the following code. You will see that all columns besides the first four will contain NA values. These are the values that must be manually input to complete the *.csv*. Do not keep any NA values in your final *.csv*. Rather, erase NA's for cells that will be left blank. The following code will save the gating template created with `templateGen()` as *gt.csv* to your computer. This will save to your current working directory. If you are unsure of this location, type `getwd()` into the R Console to retrieve the path.

```
write.csv(gt, "gt.csv")
```

If you choose to create the gating template manually, the same conventions must be followed. Start with a blank spreadsheet. Next, fill in the 10 required column names. From there, you can use the manual gating hierarchy plotted in the previous chapter to fill in each cell population *alias*. Fill in the remainder accordingly.

There will likely be troubleshooting involved in this process. This is a great place to start if you're seeking more information on the gating template. The [openCyto GitHub](#) page is also very responsive to issues posted.



## Chapter 5

# Automate Gating

### 5.1 Load .csv into R

As noted in the previous chapter, there is a sample gating template titled “*partial.csv*” with the sample data. This may serve as a guide to creating your own. When the .csv gating template is complete, it is then read into R and saved as `gt`, a `gatingTemplate` object.

```
gt <- gatingTemplate("../tutorial/partial_gt.csv")
```

### 5.2 Read in raw FCS files

Now that the `GatingTemplate` object has been loaded into R, you will need to also load raw FCS files to perform the automated gating on. For gating, these files must be in a `GatingSet` object type. The path is saved as a character matrix of file names using `list.files()`. Here, this matrix is saved as `fcs_files`. Next, `read.ncdfFlowSet()` will save FCS files as a `ncdfFlowSet` object, here named `ncfs`. The `GatingSet` function will then save the FCS files as a `GatingSet` object, here named `gs_auto`. In this form, the FCS files can be input and gated. Use the following code to perform these steps.

```
fcs_files <- list.files(path = "../tutorial/group1_v_group2", full.names = TRUE)
ncfs <- read.ncdfFlowSet(files = fcs_files)
gs_auto <- GatingSet(ncfs)
```

### 5.3 Apply Gating


At this point, you now have `GatingTemplate` and `GatingSet` object to be used for gating. Apply your `GatingTemplate` object to the `GatingSet` object, where `x` = `GatingTemplate` object and `y` = `GatingSet` object. Here, the `GatingTemplate` object is named `gt` and the `GatingSet` object is named `gs_auto`.

```
gating(x = gt, y = gs_auto)
```


### 5.4 Plot Automated Gating

Just as before, plot both the gating hierarchy and the automated gates.

```
plotGate(gs_auto[[1]])
```



```
plot(gs_auto[[1]])
```



## 5.5 Population Statistics

Both counts and frequencies of cell types in samples can be generated for analysis. Statistics can be generated from R or pulled directly from flowJo. To pull from flowJo, add `flowJo=TRUE` to either code chunk below.

### Counts

```
head(getPopStats(gs_auto, statistic="count"))
```

```
##           name Population      Parent Count ParentCount
## 1: X_group1_1 nonDebris      root  87715      142158
## 2: X_group1_1 singlets nonDebris 80775       87715
## 3: X_group1_1 lymph singlets 69718       80775
## 4: X_group1_1 viable lymph 40986       69718
## 5: X_group1_1 CD3 viable 40941       40986
## 6: X_group1_1 CD8 CD3 1888       40941
```

### Frequencies

```
head(getPopStats(gs_auto, statistic="freq"))
```

```
##           name Population      Parent Count ParentCount
## 1: X_group1_1 nonDebris      root  87715      142158
## 2: X_group1_1 singlets nonDebris 80775       87715
## 3: X_group1_1 lymph singlets 69718       80775
## 4: X_group1_1 viable lymph 40986       69718
## 5: X_group1_1 CD3 viable 40941       40986
## 6: X_group1_1 CD8 CD3 1888       40941
```



## Chapter 6

# Plot Customization

Below are two common customizations that you may want to include when plotting gates.

1. Adjust plot axes
2. Binning in a plot

### 6.1 Adjusting plots


#### 6.1.1 Adjust plot axes

As seen in Chapter 2, it may be necessary to adjust the plot axes in order to best view the gates. This is done using the code below. Setting `xlim` and `ylim` to “data” adjusts plot based on the actual data range, rather than instrument specifications. Custom ranges can also be input numerically.


```
flowWorkspace.par.set("plotGate", list(xlim = "data",  
                                       ylim = "data"))
```

Here is a comparison of `xlim` and `ylim` set as “instrument” and then “data”.

**Instrument**



## Data




### 6.1.2 Binning

Adjusting the bins changes will alter the resolution of each plot. The default is `xbin = 32`. Below


#### Default Binning

```
plotGate(gh)
```



### Greater Resolution

```
plotGate(gh, xbin = 128)
```





## Chapter 7

# Appendix

### 7.1 Function Definitions

| Function Name                 | Use  |
|-------------------------------|--|
| <code>openWorkspace()</code>  | function used to read in wsfile                                      |
| <code>parseWorkspace()</code> | function to extract FCS files from ws                                |
| <code>templateGen()</code>    | function used to generate a .csv template from existing manual gates |
| <code>gatingTemplate()</code> | function used to read in .csv template                               |
| <code>gating()</code>         | function used to apply gates to a gating set                         |
| <code>plot()</code>           | function to visualize gating tree                                    |
| <code>plotGate()</code>       | function to visualize gates  |

### 7.2 R Object Definitions

| R Object Name           | Use                                     |
|-------------------------|---|
| <code>wsfile</code>     | flowJo .wsp file location               |
| <code>ws</code>         | read in data from flowJo                |
| <code>gating_set</code> | parsed FCS files to be gated            |
| <code>gh</code>         | subset of <code>gating_set</code>       |
| <code>gt</code>         | .csv gating template                    |
| <code>gs_auto</code>    | raw FCS files to apply automated gating |

### 7.3 CSU Flow Cytometer Information

The flow cytometry equipment at CSU does not require further compensation of data. Other tutorials may highlight the steps to compensate data, but these steps are not necessary to CSU at this moment. In the event that equipment changes, it may be necessary to complete compensation steps to prepare data. More on the current equipment used as CSU here.

### 7.4 Transform data for better visualization

Transformation of these data may allow for better plotting. One common form of transformation for flow cytometry analysis is biexponential.



## Chapter 8

# References

1. Data File Standard for Flow Cytometry. (n.d.). International Society for Advancement of Cytometry. Retrieved May 10, 2019, from [http://software.broadinstitute.org/cancer/software/genepattern/attachments/fcs\\_3\\_1\\_standard.pdf](http://software.broadinstitute.org/cancer/software/genepattern/attachments/fcs_3_1_standard.pdf)
2. Duggan, R. (2014, August 20). Cytometry on Air: Analyzing Flow Cytometry Data using R. [https://www.youtube.com/watch?v=\\_B7mo6dB3BU](https://www.youtube.com/watch?v=_B7mo6dB3BU)
3. Finak, G. (2014, July 10). OpenCyto Practical Workshop. Retrieved February 14, 2019, from <https://www.bioconductor.org/help/course-materials/2014/BioC2014/OpenCytoPracticalComponent.html>
- 4 How to plot gated data. (n.d.). Retrieved March 10, 2019, from <https://www.bioconductor.org/packages/devel/bioc/vignettes/flowWorkspace/inst/doc/plotGate.html>
5. How to write a csv gating template. (n.d.). Retrieved March 3, 2019, from [https://www.bioconductor.org/packages/devel/bioc/vignettes/openCyto/inst/doc/HowToWriteCSVTemplate.html#14\\_gating\\_method\\_that\\_generates\\_multiple\\_populations](https://www.bioconductor.org/packages/devel/bioc/vignettes/openCyto/inst/doc/HowToWriteCSVTemplate.html#14_gating_method_that_generates_multiple_populations)
6. Saving your Analysis. (n.d.). Retrieved May 5, 2019, from <http://docs.flowjo.com/vx/workspaces-and-samples/ws-savinganalysis/>
7. Verschoor CP, Lelic A, Bramson JL and Bowdish DME (2015) An introduction to automated flow cytometry gating tools and their implementation. *Front. Immunol.* 6:380. doi: 10.3389/fimmu.2015.00380